Quantitative Systems Pharmacology Modeling of Loncastuximab **Tesirine Combined With** Mosunetuzumab and **Glofitamab Guides Dosing in B-Cell Lymphoma**

Yuezhe Li,¹ A. Katharina Wilkins,¹ Timothy Knab,¹ Joseph P. Boni²

¹Metrum Research Group, Tariffville, CT; ²ADC Therapeutics America, Murray Hill, NJ Presenter: Matthew Riggs, Metrum Research Group, Tariffville, CT, USA





pies of this poster obtained through the uick Response (QR) Code are for personal use only and may not be reproduced without permission from the author of this poster.

Poster presented at the 2024 Annual Meeting of the American College of Clinical Pharmacology. September 8-10, 2024, Bethesda, MD, USA | Poster originally presented at the American Association for Cancer Research Annual Meeting 2024. April 5-10, 2024, San Diego, CA, USA

CONCLUSIONS

- Substantially more tumor growth inhibition (TGI) was predicted to occur at the end of cycle 3 (C3) of loncastuximab tesirine (loncastuximab tesirine-lpyl [Lonca]) + mosunetuzumab (Mosun) combination therapy and at the end of C4 of Lonca + glofitamab (Glofit) than with any of these agents alone
- Increased doses of Lonca in combination therapy regimens were predicted to have limited additional benefit to TGI; however, additional cycles of combination therapy were predicted to enhance TGI
- Simulations suggest that Lonca doses could be reduced to allow for improved tolerability for longer periods of time, up to the point of maximal benefit, although clinical testing is needed to explore these findings
- The ongoing LOTIS-7 clinical study (NCT04970901) evaluates the safety and anticancer activity of Lonca + Mosun and Lonca + Glofit

INTRODUCTION

- CD19 and CD20, B-lymphocyte surface antigens, are clinically validated therapeutic targets for the treatment of B-cell malignancies^{1,2}
- Lonca is an antibody-drug conjugate comprising an anti-CD19 antibody conjugated to a pyrrolobenzodiazepine dimer cytotoxin that is approved as monotherapy for heavily pretreated patients with relapsed/refractory (R/R) diffuse large B-cell lymphoma (DLBCL)³
- Mosun and Glofit are CD20 × CD3 T-cell–engaging bispecific antibodies (bsAbs) that redirect T cells to eliminate malignant B cells^{4,5} - Mosun and Glofit both target a different B-cell surface antigen than Lonca; hence, combining Lonca with either of these bsAbs is expected to result in efficacy beyond that of monotherapy
- Previously, a novel physiologically based pharmacokinetic (PBPK)–quantitative systems pharmacology (QSP) model was developed⁶ and validated with Lonca monotherapy clinical observations in patients with R/R DLBCL⁷

OBJECTIVE

• To ascertain optimal Lonca + Mosun and Lonca + Glofit administration and to understand important determinants of exposure leading to the optimal dosing regimen in the context of potential CD19 expression heterogeneity

METHODS

Model Construction

• The Lonca QSP model was based on the previously validated and published Lonca PBPK-QSP model.^{6,7} This model was reduced in physiological complexity to be compatible with a published QSP model of Mosun⁸ while maintaining the core functionality of predicting tumor dynamics during Lonca monotherapy (**Figure 1**)

Figure 1. (A) Integration of Mosun subcutaneous administration into the PK-QSP model and (B) depiction of the QSP model for Lonca + Mosun/Glofit combination therapy



The following parameters were used to model Mosun subcutaneous administration: bioavailability = 0.7; absorption rate = 0.26 day-1. Glofit, glofitamab; IL-6, interleukin 6; Lonca, loncastuximab tesirine; Mosun, mosunetuzumab; PK, pharmacokinetics; QSP, quantitative systems pharmacology.

- The following assumptions were made for integrating the PKPD models:
- Lonca or Mosun/Glofit can induce tumor B-cell killing, though Lonca-killed tumor B cells would linger in the "dying" state. The rate for removing dying tumor B cells was modulated to account for different surface CD19 copy numbers
- The tumor is composed of T cells and malignant B cells (i.e., no healthy B cells)
- Tumor volume is calculated based on the tumor B-cell number (the total of tumor B cells and dying tumor cells)
- T cells can enter or leave the tumor as they enter or leave other tissues

TGI Modeling for Virtual Population

- The following assumptions and ranges were used to describe potential responses to Lonca + Mosun/Glofit combination therapy:
- ≥80% of the tumor cells are CD19⁺/CD20⁺; the remaining tumor cells are CD19^{-/low}/CD20⁺
- CD19^{-/low}/CD20⁺ tumor cells represented cells with 0-1 CD19 epitopes per cell, resulting in varied susceptibility to Lonca from completely insensitive to somewhat sensitive
- The Lonca-induced maximum killing rate of CD19+/CD20+ tumor cells varied by ±10%
- Initial tumor volume = 1-10 cc
- Malignant B-cell proliferation rate = 0.0-0.15 day⁻¹
- Lonca was given every 3 weeks (Q3W) at the following doses:
- 150 µg/kg for 2 doses followed by 75 µg/kg
- 120 µg/kg for 2 doses followed by 75 µg/kg
- 90 µg/kg
- Mosun 5 mg on cycle 1 day 1 (C1D1), 45 mg on C1D8, 45 mg on C1D15, and 45 mg Q3W for C2+
- Glofit 2.5 mg on C1D8 following obinutuzumab pretreatment on C1D1, Glofit 10 mg on C1D15, and Glofit 30 mg Q3W for C2+

TGI Modeling Parameters for a Prototypical Patient

- Initial tumor volume = 49.5 cc
- The tumor was composed of 100% CD19⁺/CD20⁺ malignant B cells
- Initial T-cell:B-cell ratio in tumor = 0.06
- Tumor B-cell proliferation rate = 0.03 day⁻¹ (tumor doubling time \approx 23 days)
- Tissue B cells and T cells were pre-equilibrated, and no depletion was assumed

RESULTS

QSP Modeling of Lonca + bsAbs

• Virtual population pharmacodynamic modeling results for Lonca + Mosun (Figure 2) and Lonca + Glofit (Figure 3) both indicated that the addition of Lonca outperformed and showed greater depth of response than either of these bsAbs alone - Both simulations indicated that TGI potential was independent of Lonca dosage

Figure 2. Virtual population simulations of tumor responses to Lonca + Mosun versus Mosun monotherapy



The solid line represents the median of the simulation, and the shaded area represents the 95% CI. CR, complete response; Lonca, loncastuximab tesirine; Mosun, mosunetuzumab; PD, progressive disease; PR, partial response; Q3W, every 3 weeks; SD, stable disease.

• Tumor volume reduction with Mosun monotherapy was predicted to plateau after day 28, whereas Lonca + Mosun combination therapy was predicted to continue to provide TGI (**Figure 4**). Similar results were observed for Glofit monotherapy compared with Lonca + Glofit combination therapy (data not shown)

Figure 4. Long-term tumor dynamics in response to Lonca + Mosun versus Lonca or Mosun monotherapy in a prototypical patient receiving both at their prescribed doses



Prototypical patient received Mosun doses of 5 mg on C1D1, 45 mg on C1D8 and C1D15, followed by 45 mg Q3W, and/ or Lonca doses of 150 µg/kg for the first 2 cycles, followed by 75 µg/kg Q3W. This patient's tumor was homogeneous (100% CD19+/CD20+). C, cycle; D, day; Lonca, loncastuximab tesirine; Mosun, mosunetuzumab; Q3W, every 3 weeks.

Presenter information

Matthew Riggs: mattr@metrumrg.com

Acknowledgments

The analysis was funded by ADC Therapeutics SA and partially funded by Sobi. Medical writing and editorial support, provided by Citrus Health Group, Inc. (Chicago, IL), was provided in accordance with Good Publication Practices (GPP 2022) and funded by ADC Therapeutics SA and Sobi.

Disclosures

Y Li: employee of Metrum Research Group. AK Wilkins: employee of Metrum Research Group. T Knab: employee of Metrum Research Group. JP Boni: was an employee of ADCT Therapeutics at the time of the study and is a current equity holder at ADC Therapeutics SA. **M Riggs (Presenter)**: employee of Metrum Research Group.

References

- 1. Bailly S, et al. *Hematol Oncol*. 2022;40(4):505-517. 2. Salles G, et al. Adv Ther. 2017;34(10):2232-2273.
- 4. Budde LE, et al. Lancet Oncol. 2022;23(8):1055-1065.
- 5. Dickinson MJ, et al. N Engl J Med. 2022;387(24):2220-2231.
- 6. Utsey K, et al. Clin Pharmacol Drug Dev. 2023;12(12):123-125.
- 7. Caimi PF, et al. *eJHaem*. 2023:5(1):76-83.
- 8. Hosseini I, et al. NPJ Syst Biol Appl. 2020;6(1):28.

Glofit monotherapy PR/CR Lonca dose (Q3W) =___ 150/75 μg/kg 💻 120/75 μg/kg Lonca + Glofit 💻 90 μg/kg — 0 μg/kg 0 7 14 21 28 35 42 49 56 63

Figure 3. Virtual population simulations of tumor responses to Lonca + Glofit versus Glofit monotherapy

The solid line represents the median of the simulation, and the shaded area represents the 95% CI. CR, complete response; Glofit, glofitamab; Lonca, loncastuximab tesirine; PD, progressive disease; PR, partial response; Q3W, every 3 weeks; SD, stable disease.

Time, day

• When examined by the total dosage of Lonca received, the TGI provided by Lonca + bsAb combination therapy was predicted to increase with additional cycles of Lonca while being independent of the Lonca dosage level (**Figure 5**)





The horizontal black dashed line represents the no response/partial response separation. Simulated dose of Glofit: 2.5/10/30 mg. Glofit, glofitamab; Lonca, loncastuximab tesirine.